FLUORESZENZLICHT- ENERGIE

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INHALT

Was ist Licht?

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
WAS IST LICHT

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
The light emitted and the fluorescence generated have different benefits on the treated skin

SHORT PENETRATION


- Control of *P. acnes* bacterial colonization
- Reduction of microbial induced inflammation
- Collagen stimulation

DERMAL TISSUE PENETRATION

- Fibroblast activation & proliferation
- Stimulate healing response
- Skin rejuvenation

DERMIS TISSUE PENETRATION

- Vascular activation
- General inflammatory reduction
- Collagen induction and rejuvenation

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The background of the slide is a photograph of three petri dishes containing yellow agar with bacterial colonies. They are placed on a blue light box, which is a common piece of laboratory equipment used for culturing light-sensitive bacteria. The scene is dimly lit, with the primary light source being the blue glow of the light box.

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Sharing the knowledge

LITERATURANGABEN

1. Kleresca® Acne Treatment Instruction for Use
2. Kleresca® Skin Rejuvenation Instruction for Use
3. Kleresca® Rosacea Treatment Instruction for Use
4. Kleresca® Pre-Post Treatment Instruction for Use
5. Jalili, A. 2018. Chromophore gel-assisted phototherapy. A novel and promising photobiomodulation therapy for facial inflammatory skin diseases and skin aging. J Ästhet Chir. 2018. Early online 20 February
6. Nikolis, A. *et al.* An extension of a multicenter, randomized, split-face clinical trial evaluating the efficacy and safety of chromophore gel-assisted blue light phototherapy for the treatment of acne. Int. J. Dermatol. 2017; Accetped
7. Antoniou, C.*et al.* . A multicenter, randomized, split-face clinical trial evaluating the efficacy and safety of chromophore gel-assisted blue light phototherapy for the treatment of acne. 2016; Int J Dermatol, 2016, 55; 1321 – 1328
8. Nikolis, A., *et al.* . A randomized, placebo-controlled, single-blinded, split-faced clinical trial evaluating the efficacy and safety of KLOX-001 gel formulation with KLOX light-emitting diode light on facial rejuvenation. Clinical, Cosmetic and Investigational Dermatology, 2016; 9, 115–25.
9. Scarcella *et al.* 2018,. Treatment of solar lentigines using a combination of PICO-Laser and Biophotonic treatment. Int J Dermatol (In press)
10. Braun, S.A.& Gerber, P. A photoconverter-gel assisted blue light therapy for the treatment of rosacea. Int. J Dermatol. 2017 [Epub ahead of print]
11. Sannino, M. *et al.* A novel platform for treating rosacea subtypes 1, 2 and 3 using fluorescent light. Int J Dermatol 2018. In press
12. Nielsen, Schoedt & Bak-Christensen. Clinical Evaluation – Kleresca® Biophotonic Treatments. 2017. Kleresca® Data on File
13. Nielsen, M.E., et al (2017). Introducing: photobiomodulation by low energy chromophore-induced fluorescent light. Mechanisms of Photobiomodulation Therapy IV, SPIE Photonics West BIOS, San Francisco, 2017; 28 January – 2 February
14. Barolet, D. Light-emitting Diodes: A Brief Review and Clinical Experience. J. Clin. Aesthet. Dermatol. 2008; 8, 36 – 44.
15. Freitas, L. F. de & Hamblin, M. R. Proposed Mechanisms of Photobiomodulation or Low-Level Light Therapy. IEEE J Sel Top Quantum Electron 22, 2016.
16. Wang, Y. *et al.* Photobiomodulation (blue and green light) encourages osteoblastic differentiation of human adipose-derived stem cells: role of intracellular calcium and light-gated ion channels. Sci. Rep 2016; 6, 33719