



INTRODUCTION TO ROSACEA

More than just redness



CONTENTS

What is rosacea?

Treating rosacea with Kleresca®
Rosacea Treatment

Take home messages



WHAT IS ROSACEA?

Rosacea: more than just redness

WHAT IS ROSACEA^{10,11}

- Chronic inflammatory skin disorder
- It affects primarily the central areas of the face (cheeks, nose, chin and central forehead)
- Often appears with periods of exacerbations and remissions of cutaneous signs



THE SIGNS OF ROSACEA^{10,11}

PRIMARY FEATURES

- Flushing (transient erythema)
- Persistent redness (non-transient erythema)
- Papules and pustules, nodules sometimes
- Telangiectasia (visible dilated blood vessels), common but not always present, more visible if redness fades

SECONDARY FEATURES

- Burning or stinging
- Plaques, larger red areas
- Dry appearance
- Oedema
- Ocular manifestations like tearing and redness
- Peripheral location, present in neck, chest, scalp, ears and back
- Phymatous changes

CLASSIFICATION OF ROSACEA^{10,11}

The standard classification system establishes the following subtypes of rosacea:



Subtype 1: Erythematotelangiectatic rosacea
Flushing
Persistent central facial erythema
Telangiectasias (visible blood vessels) common



Subtype 2: Papulopustular rosacea
Persistent central facial erythema
Swelling
Transient papules and/or pustules
Burning and stinging



Subtype 3: Phymatous rosacea
Thickening skin and enlargement
Irregular surface nodularities
Commonly rhinophyma (thickening in nose area), but can appear elsewhere



Subtype 4: Ocular rosacea
Eye irritation
Burning sensation and blood-shot eyes
Gritty feeling



CLASSIFICATION OF ROSACEA^{10,11}



Subtype 1



Subtype 2



Subtype 3



Subtype 4

CLASSIFICATION OF ROSACEA

- The different subtypes of rosacea, can have mild to severe symptoms ^{10,11}
- The latest research recommends that to be diagnosed with rosacea a patient must present²³:
 - Persistent erythema or phymatous changes
- Additionally, two or more primary features are also considered diagnostic²³
- In most cases a combination of primary and secondary features are seen in a patient, but seldomly appear all together^{10,11}



WHAT IS THE CAUSE OF ROSACEA?^{11,12}

The exact causes of rosacea are not completely understood. Latest findings suggest the following:

- Abnormal functioning of the immune system (induces inflammatory cascades)
- Neurovascular dysregulation (leads to vasodilation and development of abnormal superficial blood vessels)



LIVING WITH ROSACEA^{12,13}

- May significantly affect patient's lives, leading to psychological and social distress over appearance
- Negative effect on self-esteem and self-confidence
- Overall reduction in the patient's quality of life



TRIGGERS OF ROSACEA⁷

Rosacea symptoms can be a response to certain triggers such as:

- Anxiety and stress
- Food (especially spicy food)
- Alcohol
- Drugs
- Wind

Patients are often educated to avoid these triggers to keep their rosacea under control

ROSACEA TRIGGERS

Many patients notice certain triggers seem to make their symptoms worse.
An NRS survey of 1,066 rosacea patients found that the most common factors included ²⁴:

81%	SUN EXPOSURE	41%	INDOOR HEAT
79%	EMOTIONAL STRESS	41%	CERTAIN SKIN-CARE PRODUCTS
75%	HOT WEATHER	36%	HEATED BEVERAGES
57%	WIND	27%	CERTAIN COSMETICS
56%	HEAVY EXERCISE	15%	MEDICATIONS
52%	ALCOHOL CONSUMPTION	15%	MEDICAL CONDITIONS
51%	HOT BATHS	13%	CERTAIN FRUITS
46%	COLD WEATHER	10%	MARINATED MEATS
45%	SPICY FOODS	9%	CERTAIN VEGETABLES
44%	HUMIDITY	8%	DAIRY PRODUCTS


FACTS ABOUT ROSACEA^{11,14,18}

- Appears to be quite common
 - Prevalence rates range from **0,09%** to **22%** (depends on the study)
 - Worldwide, **5.46 %** (414,960,000 sufferers)
- Most frequently observed in patients with fair skin^{73,74}, often called the ‘Curse of the Celts’ (prevalent among Irish)⁷⁵
- Three times more common in women
- May occur at any age, typically between 30 – 50
- Rhinophyma typically appears between 50 - 70


FACTS IN IMAGES^{11,14,18}



Approximately 415 million people worldwide have rosacea



It is most frequently observed in patients with fair skin



It is 3 times more common in women than in men



It may occur at any age, typically between the age of 30 to 50



INCREASING INCIDENCE RATE¹⁵

- The steady growth and aging of the population have raised the incidence of rosacea over the past decade
- Rosacea may occur at any age, but tends to first appear during middle age
- It now affects the massive **Baby Boom generation**




COMMON BUT UNDERDIAGNOSED CONDITION^{15,16}

- The earliest signs of rosacea are often overlooked because people assume they are temporary and they will go away⁶⁹
- **95%** of rosacea patients had **known little or nothing about its signs and symptoms prior to their diagnosis**
- A study on the prevalence of rosacea in Germany and Russia found⁷²:
 - Over 80% of rosacea subjects were not previously diagnosed
 - Within the previous year, 47.5% of rosacea subjects had received no rosacea care
 - **23.7%** of rosacea subjects **had received topical and/or systemic drugs**



TREATING ROSACEA



TREATMENT OF ROSACEA¹³⁻¹⁵

- Most of the current treatments have different associated side effects (scarring, hyper/hypo pigmentation, risk of bleeding, pain...)
- They are moderately successful in controlling symptoms, but relapse usually occurs shortly after the cessation of the therapy

IMPORTANCE OF EARLY DIAGNOSIS⁷

- **Early diagnosis** and **treatment** is necessary to prevent symptoms developing from mild to moderate to severe
- Rosacea often increases in **severity** over time, especially if it is left untreated



TREATING ROSACEA WITH FLUORESCENT LIGHT ENERGY²⁻⁸

- Non-invasive biophotonic treatment
- Capacity to penetrate to various depths of the skin and to stimulate the skin tissues and cells
- Designed to target the signs and symptoms of rosacea




INSPIRED BY
PHOTOSYNTHESIS



WORKS AT
CELLULAR LEVEL



REDUCES
INFLAMMATION



FEEL THE CHANGE^{2-4,6-8}

- Reduces inflammation, reducing the presence of papules and pustules
- Reduces erythema and blushing by improving microvascularization
- Reduces the overall stress level of the skin, thereby reducing the feeling of burning and stinging
- Induces a healing response, improving the overall skin quality



KLERESCA® ROSACEA TREATMENT

- Non-invasive and comfortable treatment²⁻⁶
- High safety and efficacy²⁻⁶

Kleresca®
Photoconverter Gel



Kleresca® Rosacea Treatment

Kleresca®
Multi-LED Lamp



Rejuvenation mode

Kleresca®
Rosacea Treatment



9 min



HOW THE TREATMENT IS PERFORMED*¹



The skin is cleaned and Kleresca[®] gel is applied



The gel is illuminated for nine minutes using the multi-LED Kleresca[®] lamp, together creating fluorescent light energy that stimulates the skin



The gel is removed and the skin is cleaned and moisturised. With little to no-downtime, even make up can be applied immediately afterwards

*In general, Kleresca[®] is seen as a pleasant, comfortable and well-tolerated treatment


PLEASANT AND COMFORTABLE^{2-7,9}

- Kleresca® Rosacea Treatment has a high safety profile and is a well-tolerated treatment
- In general, the treatment is considered to be a pleasant and comfortable experience
- Side effects may include redness, swelling and hyperpigmentation (bronzing of an area of the skin)
- All side effects are transient and do not require the intervention of the clinician



REDUCTION OF ROSACEA SYMPTOMS⁷

The graphic shows a clear reduction of redness and inflammation in patients suffering with rosacea type II treated with Kleresca[®]



BEFORE AND AFTER⁶



Baseline




5 weeks after

Marked reduction of
inflammation



TAKE HOME MESSAGES



TAKE HOME MESSAGES - ROSACEA

- Chronic inflammatory skin disorder, affecting primarily central areas of the face. Often with periods of exacerbations and remissions
- **Main signs:** flushing, persistent redness, papules and pustules, telangiectasies, phymatous changes, etc. In most cases patients present a combination of primary and secondary signs
- The standard classification system establishes 4 subtypes of rosacea
 - **Subtype I:** Erythematotelangiectatic rosacea
 - **Subtype II:** Papulopustular rosacea
 - **Subtype III:** Phymatous rosacea
 - **Subtype IV:** Ocular rosacea
- Quite common, affecting almost **415 Million** sufferers worldwide, frequently in fair skin. The exact causes are not well understood
- Common but underdiagnosed, **95% of patients** had known little or nothing about its signs and symptoms before diagnosis
- May **significantly affect patients' lives**, leading to psychological and social distress over appearance
- **Early diagnosis** and **treatment** is necessary to prevent symptoms developing from mild to moderate to severe
- Can be treated with **fluorescent light energy** (Kleresca®). A non-invasive and gentle treatment designed to target its signs and symptoms
- Kleresca® Rosacea Treatment reduces inflammation, reducing the presence of papules and pustules and it reduces erythema and blushing by improving microvascularization
- It also reduces the overall stress level of the skin, reducing the feeling of burning and stinging and it induces a healing response, improving the overall skin quality





Do you have further questions?
Send us a mail at academy@kleresca.com

Do you want to check your knowledge about rosacea? Do the test to see how much you have learned and get your certificate of completion!

TAKE THE TEST



Sharing the knowledge

REFERENCES

1. Kleresca® Rosacea Treatment Instruction for Use
2. Jalili, A. 2018. Chromophore gel-assisted phototherapy. A novel and promising photobiomodulation therapy for facial inflammatory skin diseases and skin aging. *J Ästhet Chir.* 2018. Early online 20 February
3. Nikolis, A. et al. An extension of a multicenter, randomized, split-face clinical trial evaluating the efficacy and safety of chromophore gel-assisted blue light phototherapy for the treatment of acne. *Int. J. Dermatol.* 2017; Accepted
4. Antoniou, C. et al. A multicenter, randomized, split-face clinical trial evaluating the efficacy and safety of chromophore gel-assisted blue light phototherapy for the treatment of acne. 2016; *Int J Dermatol*, 2016, 55; 1321 – 1328
5. Nikolis, A., et al. A randomized, placebo-controlled, single-blinded, split-faced clinical trial evaluating the efficacy and safety of KLOX-001 gel formulation with KLOX light-emitting diode light on facial rejuvenation. *Clinical, Cosmetic and Investigational Dermatology*, 2016; 9, 115–25.
6. Braun, S.A. & Gerber, P. A photoconverter-gel assisted blue light therapy for the treatment of rosacea. *Int. J. Dermatol.* 2017 [Epub ahead of print]
7. Sannino, M. et al. A novel platform for treating rosacea subtypes 1, 2 and 3 using fluorescent light. *Int J Dermatol* 2018. In press
8. Nielsen, Schoedt & Bak-Christensen. Clinical Evaluation – Kleresca® Biophotonic Treatments. 2017. Kleresca® Data on File
9. Nielsen, M.E., et al (2017). Introducing: photobiomodulation by low energy chromophore-induced fluorescent light. *Mechanisms of Photobiomodulation Therapy IV*, SPIE Photonics West BIOS, San Francisco, 2017; 28 January – 2 February
10. Wilkin, J., et al. Standard grading system for rosacea: report of the National Rosacea Society Expert Committee on the classification and staging of rosacea. *J Am Acad Dermatol.* 2004 Jun;50(6):907-12.
11. Wilkin, J., et al. Standard grading system for rosacea: report of the National Rosacea Society Expert Committee on the classification and staging of rosacea. *J Am Acad Dermatol.* 2002 Apr;46(4):584-7.
12. Zeichner, A., et al. Quality of Life in individuals with erythematotelangiectatic and papulopustular rosacea: findings from a web-based survey. *J Clin Aesthet Dermatol.* 2018; 11(2):47 – 52
13. Say, E.M., et al. Treatment Outcomes of Long-Pulsed Nd: YAG Laser for Two Different Subtypes of Rosacea. *J Clin Aesthet Dermatol.* 2015 Sep;8(9):16-20.
14. J. Spoendlin, J. J. Voegel, S. S. Jick, C. R. Meier. (2012, September). A study on the epidemiology of rosacea in the U.K. *British Journal of Dermatology*, 167(3), 598-605.
15. National Rosacea Society. (2018). Rosacea Now Estimated to Affect at Least 16 Million Americans. Retrieved January 31, 2018, from www.rosacea.org: https://www.rosacea.org/rr/2010/winter/article_1.php
16. Mandal, A. (2013, September 22). Rosacea Epidemiology. Retrieved January 31, 2018, from www.news-medical.net: <https://www.news-medical.net/health/Rosacea-Epidemiology.aspx>
17. National Rosacea Society. (2017). Rosacea Awareness Month. Retrieved February 1, 2018, from www.rosacea.org: <https://www.rosacea.org/patients/ram>
18. Gether, L. et al. Incidence and Prevalence of Rosacea: a systematic review and meta-analysis. *Br. J. Dermatol* 2018 [Epub ahead of print]
19. J. Tan, H. Schöfer, E. Araviiskaia, F. Audibert, N. Kerrouche, M. Berg on behalf of the RISE studygroup. (2016, March). Prevalence of rosacea in the general population of Germany and Russia - The RISE study. *Journal of the European Academy of Dermatology and Venereology*, 30(3), 428-34.
20. Stanisław Jarmuda, Niamh O'Reilly, Ryszard Z'aba, Oliwia Jakubowicz, Andrzej Szkaradkiewicz and Kevin Kavanagh. (2012). Potential role of Demodex mites and bacteria in the induction of rosacea. *Journal of Medical Microbiology*, 61, 1504–1510.
21. Thomas Jansen and Gerd Plewig. (1997). Rosacea: classification and treatment. *Journal of the Royal Society of Medicine*, 144-150.
22. National Rosacea Society. Widespread Facial Disorder May Be Linked to Genetics. www.rosacea.org/press/archive/20080602.php. Retrieved May 2018.
23. Tan et al. Updating the diagnosis, classification and assessment of rosacea: recommendations from the global ROSacea COnsensus (ROSCO) panel. *Br J Dermatol.* 2017 Feb;176(2):431-438
24. National Rosacea Society. Rosacea Triggers Survey. Retrieved 25 May 2018 from: <https://www.rosacea.org/patients/materials/triggersgraph.php>

