# L'ÉNERGIE FLUORESCENTE

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# FLUORESCENT LIGHT ENERGY<sup>5-13</sup>

- Normalise cellular activity
- Improve overall skin quality
- Limited to no downtime
- Pleasant & comfortable treatment

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




Photoconverter gel absorbs blue light from the lamp


Chromophores in the gel convert the blue light into fluorescent light energy<sup>13</sup>

Photobiomodulation stimulates the skin's own repair mechanisms<sup>14,15</sup>  
Unique and dynamic photonic output, acting on the epidermis and the dermis<sup>13</sup>

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
# CONCLUSIONS

# TAKE HOME MESSAGES ABOUT LIGHT

- Light is important in our daily life and it is also used in many medical treatments, like dermatology conditions and wound healing.
- Light moves in waves and it is measured in nanometers (nm). Visible light ranges from 400 to 700 nm.
- Depending on its wavelength, light penetrates differently into the skin
- Biophotonics is the interaction of light with biological systems, and modulates several aspects of the cell biology.
- Photobiomodulation helps in the reduction of inflammation, increases collagen build-up and normalises the cell activity.
- Kleresca® mode of action is inspired by photosynthesis:
  - The photoconverter gel absorbs blue light from the lamp
  - Chromophores in the gel convert the blue light into fluorescent light energy
  - Photobiomodulation stimulates the skin's own repair mechanisms
- The light emitted and the fluorescence generated have different benefits on the treated skin like control of *P. acnes*, fibroblast activation and proliferation, stimulation of healing response and vascular activation.
- Pleasant and easy to use treatment

Do you have further questions?  
Send us a mail at [academy@kleresca.com](mailto:academy@kleresca.com)

Do you want to check your knowledge about what is light and the mode of action of Kleresca®? Do the test to see how much you have learned and get your certificate of completion!

TAKE THE TEST

*Sharing the knowledge*

# RÉFÉRENCES

1. Kleresca® Acne Treatment Instruction for Use
2. Kleresca® Skin Rejuvenation Instruction for Use
3. Kleresca® Rosacea Treatment Instruction for Use
4. Kleresca® Pre-Post Treatment Instruction for Use
5. Jalili, A. 2018. Chromophore gel-assisted phototherapy. A novel and promising photobiomodulation therapy for facial inflammatory skin diseases and skin aging. J Ästhet Chir. 2018. Early online 20 February
6. Nikolis, A. *et al.* An extension of a multicenter, randomized, split-face clinical trial evaluating the efficacy and safety of chromophore gel-assisted blue light phototherapy for the treatment of acne. Int. J. Dermatol. 2017; Accetped
7. Antoniou, C.*et al.* A multicenter, randomized, split-face clinical trial evaluating the efficacy and safety of chromophore gel-assisted blue light phototherapy for the treatment of acne. 2016; Int J Dermatol, 2016, 55; 1321 – 1328
8. Nikolis, A., *et al.* A randomized, placebo-controlled, single-blinded, split-faced clinical trial evaluating the efficacy and safety of KLOX-001 gel formulation with KLOX light-emitting diode light on facial rejuvenation. Clinical, Cosmetic and Investigational Dermatology, 2016; 9, 115–25.
9. Scarcella *et al.* 2018,. Treatment of solar lentigines using a combination of PICO-Laser and Biophotonic treatment. Int J Dermatol (In press)
10. Braun, S.A.& Gerber, P. A photoconverter-gel assisted blue light therapy for the treatment of rosacea. Int. J Dermatol. 2017 [Epub ahead of print]
11. Sannino, M. *et al.* A novel platform for treating rosacea subtypes 1, 2 and 3 using fluorescent light. Int J Dermatol 2018. In press
12. Nielsen, Schoedt & Bak-Christensen. Clinical Evaluation – Kleresca® Biophotonic Treatments. 2017. Kleresca® Data on File
13. Nielsen, M.E., et al (2017). Introducing: photobiomodulation by low energy chromophore-induced fluorescent light. Mechanisms of Photobiomodulation Therapy IV, SPIE Photonics West BIOS, San Francisco, 2017; 28 January – 2 February
14. Barolet, D. Light-emitting Diodes: A Brief Review and Clinical Experience. J. Clin. Aesthet. Dermatol. 2008; 8, 36 – 44.
15. Freitas, L. F. de & Hamblin, M. R. Proposed Mechanisms of Photobiomodulation or Low-Level Light Therapy. IEEE J Sel Top Quantum Electron 22, 2016.
16. Wang, Y. *et al.* Photobiomodulation (blue and green light) encourages osteoblastic differentiation of human adipose-derived stem cells: role of intracellular calcium and light-gated ion channels. Sci. Rep 2016; 6, 33719